




Materials:

Filling tube: PVC
 Discharge tube: silicone
 Reservoir bottle: PE-LD
 Base: PE-HD
 Fitting, grey: Polyamide
 Micrometer screw: brass/PP

- Squeezing bottle allows rapid filling
- Automatic zeroing
- Micrometer screw allows fine titration
- Rough titration by releasing the micrometer screw and using the press button

Note!  Burettes not in use for longer periods should be drained and tubing pressure released by opening the stopcock.

Automatic burettes, Dr. Schilling pattern

SILBERBRAND

Burette: AR-GLAS®. Error limits according to class B, DIN EN ISO 385. Calibrated to deliver (TD, Ex). Schellbach stripe. Automatic zeroing. Pack of 1.

Capacity ml	Subdivision ml	Error limit ± ml	With bottle ml	Total height mm	Cat. No.
10*	0.05	0.05	500	530	237 53
15*	0.1	0.08	500	510	237 55
25*	0.1	0.08	1000	620	237 56
50*	0.1	0.10	1000	830	237 58
25*, amber glass	0.1	0.08	1000, amber	650	237 66
50, amber glass	0.1	0.10	1000, amber	900	237 68

* reduced distance between subdivision marks

Spare parts for automatic burettes, Dr. Schilling pattern

Spare burettes

Pack of 1.

Capacity ml	Length mm	Cat. No.
10*	340	237 13
15*	320	237 15
25*	390	237 16
50*	600	237 18
25*, amber glass	420	237 23
50, amber glass	670	237 24

* reduced distance between subdivision marks

Bases

Pack of 1.

For bottle ml	Ø bottle mm	Cat. No.
500	75	237 25
1000	94	237 28

Reservoir bottles

PE-LD. Narrow-mouth bottle with tube bushing. Pack of 1.

Capacity ml	Ø mm	Height mm	Thread GL	Cat. No.
500	75	180	25	1290 55
1000	94	212	28	1290 60
1000, amber	94	212	28	1302 60

Fittings

Micrometer screw with push-button, incl. support clamp. Pack of 1.

For capacity ml	for reservoir bottle ml	Cat. No.
10	500	237 45
15	500	237 46
25 - 50	1000	237 48

Glass discharge tip

AR-GLAS®. With silicone tubing. Pack of 10.

Clear glass

Cat. No.	115 00
----------	---------------

Amber glass

Cat. No.	115 05
----------	---------------